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Agata Piestrzyńska-Kajtoch , tel. 666081328, mail: agata.kajtoch@iz.edu.pl

Grzegorz Smolucha , tel. 666081204, mail: grzegorz.smolucha@iz.edu.pl

THE INSTRUCTION OF COLLECTING MATERIAL FOR PRNP GENOTYPING (PB09)

The blood should be collected by vet or the person authorized

The person collecting the samples is obliged to:

- sterile sample collection (in case of taking several samples from different animals it is necessary to remember about changing gloves or disinfecting hands after taking material from each next animal in order to protect samples from contamination; contamination of biological material with other biological material makes it **NOT suitable** for further analysis!!!);
- animal identity verification;
- legible and durable description of the samples taken;
- checking the compliance of the labeling of the collected material with the collection protocol;
- storing and transporting samples in accordance with the following recommendations;
- signing the protocol for the collection of biological material, confirming the collection of the material in accordance with the requirements contained in the manual or presenting its own sampling procedure on the order.

Sample collecting:

1. Blood (1 ml minimum) should be collected to hermetic test-tubes with EDTA.
2. The Test-tube with blood and EDTA should be mixed by several times hand rotating.
3. Every test-tube must have the readable description containing the sample number or animal identification number (concordant with the sample collection protocol).
4. Test-tubes can be delivered personally or by post; delivery address:

Contact person*

**Laboratory of Molecular Genetics
Department of Animal Molecular Biology
National Research Institute of Animal Production
Krakowska 1
32-083 Balice
Poland**

Along with the samples, a signed order and a sample collection protocol (originals) should be sent.

The material can be delivered to the laboratory between 7.30 - 14.00.

Lack of order or incomplete documentation makes it impossible to carry out the study.

Notice:

Test-tubes should be plastic with durable, tight cap or screw cap. Do not send glass test-tubes or test-tubes with needles. Such a tubes are dangerous for people who take care of packages. They leak very often, and because of that they may destroy other packages and shipping companies can have claims. **The laboratory can decline to receipt material in damaged test-tubes or with needles.**

The time of sample transport is important for analysis quality and correctness, so it is better to send samples by courier (not by regular post). The fresh samples should be delivered to the Lab in 24-48 hours after collection. If it is impossible, the samples should be properly protected:

- Shipment on the second or third day after collection: samples should be stored in the refrigerator (2-8°C) and shipped.
- Shipment on the third day after collection or later: samples should be frozen in the freezer and shipped by courier in the package (Styrofoam) containing ice (e.g. plastic bottle with frozen water) to avoid thawing.
- The temperature conditions have an influence on the sample quality. If the samples are sent during hot weather, they should be properly protected e.g. shipped by courier in the package (Styrofoam) containing ice (e.g. plastic bottle with frozen water).
- Test-tubes should be properly protected in the package against bumps (with e.g. bubble wrap, styrofoam mouting, styrofoam marbles, paper fritters, newspaper). They shouldn't bump to each other during transportation. It happens, that insufficiently protected test-tubes crash, fracture, burst or their lids open. In such a case, samples are lost or contaminated – that's make the analysis impossible. **In such a case the Laboratory has the right to refuse the test.**

***Contact person:**

Agata Piestrzyńska-Kajtoch , tel. 666081328, mail: agata.kajtoch@iz.edu.pl

Grzegorz Smolucha , tel. 666081204, mail: grzegorz.smolucha@iz.edu.pl